

ALS Laboratory Learning

Synthetic sample handling guide | Manual demo v1

1. Identify the sample

Use the synthetic sample code DEMO-001.

Match the specimen label to the request before processing.

If a label is missing or inconsistent, stop and request clarification.

2. Prepare safely

Wear the protective equipment required for the teaching exercise.

Keep the work surface clean and use a clearly labelled container.

This demonstration contains no patient data or clinical specimen.

3. Document the check

Record the sample code, date, and the observed label match.

Separate an observation from an interpretation.

Never invent a result to complete an incomplete worksheet.

Practice question

What should be checked before processing a sample?

Answer: its identity and label. Explain any discrepancy first.